

SA-15K INTIALIZATION PROCEDURE

- 1. Press and hold the CAL key,press ON/OFF key then release both keys...."88888's appear"*
- 2. Within 3 seconds press and hold both HOLD and ZERO keys,press CAL key then release all keys....."SETTING appears."*
- 3. Press and "hold" the CAL key until the display shows "INIT A then SETTING".*
- 4. Press and hold the HOLD key followed immediately by the TARE key.C11 0 is displayed and remains unchanged.*
- 5. Press the 0 key to move to C12 0 ,press the 00 key to change to C12 1.*
- 6. Press the ENTER key and display shows "SETTING".Press the ON/OFF key to return to normal weighing mode.Proceed with calibration.*
- 7. Press the CAL key,then press the ZERO,followed by the ENTER key,CAL 0 is displayed.Wait until stable light is shown then press the ENTER key,CAL F is displayed.Apply 15kgs to the pan and wait for the stable light,then press the ENTER key,END is displayed.*
- 8. Press the CAL key to return to normal weighing mode.*