



HP Series –Span Calibration

SPAN CALIBRATION

1. Verify that the display is on.
2. Press & hold the ZERO key. “Cal 0” displayed.
3. Ensure the pan is empty & press the ZERO key. → “CAL F” displayed.
4. Place mass A on the weigh pan and press the ZERO key. “CAL END” displayed.
5. Press ON/OFF key twice to return to normal weigh mode.

MODEL	MASS A
HP12K	10 kg
HP20/22K	20 kg
HP30K	30 kg
HP40K	40 kg